

SENTINEL-CXR


Screenshots of the deployed system

Final Group Project · Deep Learning (MAIB AI 114) · S P Jain School of Global Management, Dubai
Krishna Mathur AS25DXB018 · Atharva Soundankar AS25DXB020 · Yash Petkar AS25DXB021

Captured from <https://sentinel-cxr.vercel.app> on 12 August 2026, 20:44 UTC.

These are the live application, not mockups. The system can be opened directly;
the console's demo requires no account.

Research prototype. Not a medical device.



Scrolling assembles the thorax and reveals the measured metrics

SENTINEL • CXR

MethodDashboard

Console

0.9004
empirical coverage

151 ms
on 0.1 CPU

112,120
radiographs

83
tests passing

SCROLL TO ASSEMBLE

• Consolidation

• Pneumonia

• Mass

• Effusion

• Cardiomegaly

• Infiltration

• Nodule

• Atelectasis

• Pleural Thickening

• Fibrosis

• Emphy

MEASURED, NOT CLAIMED

0.9004

Empirical coverage

against a 0.90 nominal target

143

Milliseconds

ONNX fast path, 0.1 CPU

112,120

Radiographs

30,805 patients, patient-disjoint splits


83

Tests passing

coverage, grounding, triage, fairness

Every figure above was produced by running the system, not estimated. Classification metrics on the full dataset are marked outstanding in the report rather than filled in.

https://sentinel-cxr.vercel.app · 12 August 2026, 20:44 UTC



SENTINEL · CXR

Sign in

Studies are private to your account.

Email

Password

Sign in

OR

Explore the demo — no sign-up

Opens a private sandbox pre-loaded with five studies covering triage, abstention, progression and a rejected upload. Expires after 24 hours.

No account? Create one

The demo needs no account and is the fastest way in. Note that on free-tier hosting the database is wiped whenever the backend restarts, so registered accounts may not persist between visits — the demo reissues itself automatically.

Research prototype for MAIB AI 114. Not a medical device. Do not upload identifiable patient data.

SENTINEL - CXR

REDUCED

Demo Reviewer

Sign out

Patient reference

PT-001

Groups studies into a timeline. Use a pseudonym, never a real identifier.

Upload radiograph

Ask about your studies

show me abstained studies

Ask

abstained studies

STAT cases

worsening patients

above 0.8

WORKLIST · 4

STAT

04:00

DEMO-014

Pneumothorax · 0.94

URGENT

47:00

DEMO-031

Abstained

URGENT

22:00

DEMO-007

Effusion · 0.87

ROUTINE

03:10

DEMO-007

Effusion · 0.63


No study selected

Upload a frontal chest radiograph, or choose one from the worklist.


Findings, explainability, progression, uncertainty and the drafted report appear here once a study is open.

Research prototype for MAIB AI 114. Not a medical device. Every output requires radiologist review.


https://sentinel-cxr.vercel.app · 12 August 2026, 20:44 UTC



Measured coverage: 9 of 14 labels below the 0.90 target



The fairness audit fails its own tolerance at 0.2149



SENTINEL - CXR

Dashboard

Console

Demo Reviewer

Overview

Model

Fairness

Audit

Audit record

5 of 5 studies. Click a column to sort, a row to open it.

Search...

CSV

Patient ▾	Decision ▾	Triage ▾	Latency ▾	When ▾	Findings
DEMO-014 · 1	answered	STAT	1039ms	08:39 PM	Pneumothorax
DEMO-099 · 1	rejected	ROUTINE	1503ms	08:34 PM	—
DEMO-007 · 2	answered	URGENT	1294ms	08:21 PM	Effusion, Cardiomegaly
DEMO-031 · 1	abstained	URGENT	1276ms	07:56 PM	Consolidation, Pneumonia, Infiltration, Atelectasis, Mass, Effusion, Edema...
DEMO-007 · 1	answered	ROUTINE	745ms	05:33 PM	Effusion, Cardiomegaly

Recent activity

Newest first. Click any study to open it in the console.

DEMO-014

Pneumothorax

1039ms

DEMO-099

rejected at the gate

1503ms

DEMO-007

Effusion

1294ms

DEMO-031

abstained — routed to a human

1276ms

DEMO-007

Effusion

745ms

Research prototype for MAIB AI 114. Not a medical device. Every output requires radiologist review.